

# WORK INSTRUCTION BREAKDOWN SHEET

Operation: Scanning Electron Microscope (SEM) **Operations**

Instrument: ASPEX EXplorer SEM

IMPORTANT STEPS	KEY POINTS	REASONS WHY
A logical segment of the operation when something happens to advance the work.	Anything in a step that might: 1. Make or break the job 2. Injure the worker 3. Be a Cultural Consideration 4. Make the work easier to do (i.e., “knack”, “trick”, special timing, or bit of special information).	Reasons for each key point.


# VACUUM CONTINUITY POLICY FOR SEM (Read before proceeding!)

- **Sample Chamber** is under high vacuum.
- Computer and SEM are to be left **ON at all times** to insure continuity of vacuum per NanoLab policy.
- Rapid shift to ambient atmospheric pressure can cause significant damage to the instrument’s components!
- Vacuum is maintained in **Sample Chamber** to minimize contamination and to decrease rate of oxidation of the tungsten (W) **Filament**; thus, vacuum continuity maximizes the **Filament** lifetime.

## Imaging

Imaging consists of:


- Navigation
- Brightness & Contrast (i.e., exposure)
- Magnification
- Focus (i.e., sharpness)
- Stigmator function
- Degauss function
- Saving Images



### Navigation - Perception Suite User Documentation § 3-10

Manipulate the sample location to begin imaging.

- CTRL & LEFT-CLICK to position the **Stage** beneath the electron beam.
  - Electron beam position is approximated by black circle with white dot in center.
  - Operator can select color of circle by clicking on "**Cursor**" button.
- Both of the following methods will center the **Stage** wherever the image is clicked:
  - To move between samples on the **Sample Stub Holder**, CTRL & LEFT-CLICK on the image on the left monitor.
  - To move around an individual sample, CTRL & LEFT-CLICK on the image on the right monitor.




## Brightness & Contrast (i.e., exposure) - Perception Suite User Documentation § 3-12 – 3-15

Adjust Brightness and Contrast using the labeled knobs under the **Selector/SED/QuadBSED** tabs.

Adjustments to Brightness and Contrast knobs in the imaging modes (i.e., **Selector**, **SED**, or **QuadBSED**) must be made while the respective mode is active; adjustments made in one mode do not automatically adjust the other two modes.

- **SED** (Secondary Electron Detector) emphasizes topology or structure of the sample.
- **QuadBSED** (Quadrant BackScatter Electron Detector) emphasizes compositional information.
- **Selector** allows the Operator to blend SED and QuadBSED modes.


- Unlike focusing an image, brightness and contrast are subjective. What looks “right” may depend on the intended use of the image.
- For more information on detection modes, see sections 3-13 - 3-15 in the Perception Suite User Documentation.



## Magnification - Perception Suite User Documentation § 3-17

ZOOM the image in and out to analyze areas of interest.


- Increase or decrease magnification using either the + and – buttons, or the knob.
- Image resolution significantly decreases at magnification greater than 25,000X.




## Focus (i.e., sharpness) - Perception Suite User Documentation § 3-19

Focus/sharpen the image using the **Focus** control grid.

- LEFT-CLICK and drag the cursor to the right or left to increase or decrease the focus distance.
- Focal changes may be easier to observe on the edges of objects' features.
- Focus precision is directly proportional to magnification – higher zoom levels correspond to finer focal adjustments.
- NOTE: magnification levels will automatically change when focus is adjusted.



### Stigmator - Perception Suite User Documentation § 3-19

<p>OPTIONAL: Operator may adjust the <b>Stigmator</b> by LEFT-CLICKing and dragging the cursor to slide the grid in the x and y directions individually.</p>	<p>Only If Stigmation adjustments are desired by Operator:</p> <ul style="list-style-type: none"> <li>• Manual adjustment of <b>Stigmator</b> should be performed along one axis (i.e., x or y) at a time; dragging <b>Stigmator</b> grid diagonally does not achieve desired results.</li> <li>• A poorly stigated image may display sharp edges along one axis, but fuzzy edges along another.</li> <li>• The <b>AutoStig</b> button typically yields satisfactory results.</li> <li>• Cycling between <b>Focus</b> and the <b>Stigmator</b> may be necessary a few times to achieve optimal results.</li> </ul>	<ul style="list-style-type: none"> <li>• At magnification greater than 2500x, <b>Focus</b> alone may not be sufficient to produce a sharp image.</li> <li>• The <b>Stigmator</b> provides a mechanism to fine tune image sharpness by compensating for shape imperfections inherent in electron optics.</li> <li>• <b>Stigmator</b> adjusts the shape of the electron beam from an ellipse (undesired) to a perfect circle (desired).</li> </ul> 
--	--	--

### Degauss - Perception Suite User Documentation § 3-19

<p>OPTIONAL: Operator may LEFT-CLICK the <b>Degauss</b> button for additional image improvement.</p>	<p>Only If Degaussing is desired by Operator:</p> <ul style="list-style-type: none"> <li>• Click the <b>Degauss</b> button located in the <b>Stigmator</b> control area.</li> <li>• The image will probably go out of focus but to a lesser degree with each Degaussing.</li> <li>• Note that the image will be refocused at a shorter focal length</li> <li>• Depending on updated focus/resolution, Operator may repeat the <b>Degauss</b> process a few times (typically up to 3 times) until satisfactory focus/resolution is achieved.</li> </ul>	<ul style="list-style-type: none"> <li>• The <b>Degauss</b> function corrects discrepancies between the objective lens's current and its field strength.</li> <li>• A consistent relationship between these two variables is important to ensure reproducibility in measurements of the working distance and magnification.</li> </ul>
--	--	--




### Note on image quality

For optimal imaging results, it may be necessary to cycle between **Focus**, **Stigmator**, and **Degauss** a few times. In-depth explanation of these principles is beyond the scope of this SOP; refer to the Perception Suite User Documentation if additional information is desired.

## Saving Images - Perception Suite User Documentation § 4-28

Saving Images will entail use of the **Toolbar** located on the right monitor in addition to some functions on the **Perception Console** in the left monitor.



<p>When satisfied with an optimally focused image, increase image resolution.</p>	<p>On the left monitor:</p> <ul style="list-style-type: none"> <li>• <b>Decrease <i>Scan Speed</i></b>: LEFT-CLICK and drag the cursor to slide the vertical bar to the right end of the blue range.</li> <li>• <b>Increase <i>Image Size</i></b>: LEFT-CLICK and drag the cursor to slide the vertical bar to the right one or two notches to 1024X1024 or 2048x2048.</li> </ul>	
<p>Initiate the Save function.</p>	<p>On the right monitor:</p> <ul style="list-style-type: none"> <li>• Select the <b><i>Diskette</i></b> icon.</li> </ul> <p>On the left monitor:</p> <ul style="list-style-type: none"> <li>• It can take several seconds to minutes for the Windows “<b>Save image with context data</b>” dialog box to appear.</li> <li>• NOTE: To verify Save function is proceeding, “Save pending at end-of-frame.” will display in the <b>Message Area</b> immediately to the right of the <b>Power</b> button when the <b><i>Diskette</i></b> icon is clicked.</li> </ul>	<p>Decreasing Scan Speed and increasing Image Size improves image quality and correspondingly generates more data; therefore, computer processing time increases.</p> 
<p>Finalize the Save function.</p>	<ul style="list-style-type: none"> <li>• Click on “<b><i>Path</i></b>” in the “<b>Save image with context data</b>” dialog box and choose the appropriate file folder.</li> <li>• Name the .TIF image file and include any additional information in the “<b><i>Comments</i></b>” area.</li> </ul>	<p>All settings and “<b><i>Comments</i></b>” will be saved in the metadata of the .TIF image file and can be viewed later using Aspek Personal Image Print II or other programs.</p> 

When finished inspecting and/or saving images, proceed to the **SEM STARTUP / SHUTDOWN Work Instructions** to end your session on the SEM.